

# MRI Trigger Interface 4.4

---

Interfacing the MR Scanner with SCAN



# Table of Contents

<b>Part I MRI Trigger Interface</b>	<b>1</b>
1 Contact Information .....	2
2 Device Classification .....	3
<b>Part II Installation and Details</b>	<b>6</b>
1 Installing the MRI Trigger Interface .....	6
2 In Operation .....	7
3 Troubleshooting .....	7
<b>Index</b>	<b>0</b>

# 1 MRI Trigger Interface

## *MRI Trigger Interface User Guide*

### ***Compumedics USA, Inc.***

6605 West W.T. Harris Blvd., Suite F  
Charlotte, NC 28269

USA

Telephone: 877-717-3975 (8am – 5pm EST)



Compumedics Germany GmbH  
Heussweg 25

20255 Hamburg, Germany

Telephone: +49 40 40 18 99 41

Fax: +49 40 40 18 99 49

Internet: [sales@neuroscan.com](mailto:sales@neuroscan.com)

[techsup@neuroscan.com](mailto:techsup@neuroscan.com)

[www.neuroscan.com](http://www.neuroscan.com)



## 1.1 Contact Information

For Technical Support.....

If you have any questions or problems, please contact Technical Support through any of the following routes.

If you live outside the USA or Canada, and purchased your system through one of our international distributors, please contact the **distributor** first, especially if your system is under warranty. In all other cases, please use **techsup@neuroscan.com**, or see the other Support options on our web site (<http://www.neuroscan.com>). Or, if you live in the USA or Canada, please call **1-800 474-7875**. International callers should use **877-717-3975**.

For Sales related questions, please contact your local distributor, or contact us at **sales@neuroscan.com**.

### **Copyright © 2009 - Compumedics Neuroscan.**

All rights reserved. Printed in the United States of America. No part of this manual may be used or reproduced in any form or by any means, or stored in a database or retrieval system, without prior written permission of the company. Making copies of any part of this document for any purpose other than your own personal use is a violation of United States copyright laws. For information, contact Compumedics/Neuroscan.

3570C MRI Trigger Interface Users Guide

## 1.2 Device Classification



### **ATTENTION: CONSULT ACCOMPANYING DOCUMENTS BEFORE USING**

The MRI Trigger Interface PN3361 assembly is a 5 VOLT DC powered assembly and should be used only according to the manufacturer's instructions. Replacement parts and accessories may be obtained from the manufacturer.

**Manufacturer:** **Compumedics Neuroscan USA Ltd.**  
6605 West W.T. Harris Blvd., Suite F  
Charlotte, NC 28269 USA  
Phone: 877-717-3975 (8am – 5pm EST)

E-mail:

[techsup@neuroscan.com](mailto:techsup@neuroscan.com)  
Web: [www.neuroscan.com](http://www.neuroscan.com)

There is no guarantee that interference will not result from operation of this device in proximity or connected to some other device. If interference occurs, the user or operator is encouraged to try and correct the interference by one or more of the following measures: (1) Change the orientation of the two devices relative to one another; (2) Increase the separation between the two devices; (3) Check the power source and grounding for the two devices; and (4) Consult the dealer, Technical Support, or an experienced technician for help.

### **Classification**

The assembly is ordinary equipment not protected against ingress of water and should not be used in the presence of any spilled liquids. It is not designed to be suitable for use in the presence of a flammable anesthetic mixture of air and oxygen or nitrous oxide. The device is capable of continuous operation.

### **Technical Description**

MRI Trigger Interface Input: 5 Volts D.C. +/- .25 volts @ <30mA

#### *MRI Trigger Interface*

Weight:	75gr
Dimensions	Height: 66.5cm
	Width: 66.5cm
	Depth: 18.8cm

### **Shipping and Storage Maximum Limits**

-20°C to +70°C, 10% to 100% humidity, non-condensing RH, 500hPa to 1060hPa. After unpacking, allow devices to adjust to room temperature for at least two hours prior to interconnection and application of power.

### **Operational Limits**

+15°C to +30°C, 25% to 95% humidity, non-condensing RH, 700hPa to 1060hPa pressure.

### **Warnings and Precautions**

#### **Instructions**

Read instructions before operating the device.

The MRI Trigger Interface operates using 5 volts D.C. supplied by a host PC USB cable.

This device is not equipped with appropriate alarms required for use in monitoring clinical parameters of a patient where it is necessary to alert the user of situations which could lead to death or severe deterioration of the patient's state of health.

**CAUTION:** MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC (ElectroMagnetic Compatibility) and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.

**CAUTION:** Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT. Use of this type of equipment in close proximity to the Power Unit should be prohibited. If portable or mobile RF communications devices are used in the vicinity of the Power Unit, the user or operator should verify normal operation of the device.

### **Environment**

The MRI Trigger Interface is designed to be used in a research laboratory environment. Extremes of humidity, temperature, or pressure should be avoided. The device should not be used in a location where contact with liquids is possible, and if liquids are spilled on or in the area of the device, it should not be used until it can be ensured that the fluid or its residue will not affect device operation. Questions should be directed to the manufacturer or its representatives.

### **Cleaning Instructions**

The MRI Trigger Interface enclosure may be cleaned with a damp sponge or cloth and mild nonabrasive cleanser. Take care to ensure that liquid does not spill in or on the device. Do not use abrasives or detergents.

### **Repair**

There are no user serviceable parts in the MRI Trigger Interface. Contact your dealer or Compumedics/Neuroscan Technical Support if you believe the MRI Trigger Interface is in need of repair.

### **Maintenance**

Compumedics/Neuroscan suggests that the assembly be visually inspected at least once per year. Replace any worn or frayed cables, and contact your dealer or Compumedics/Neuroscan technical support if you have concerns about what you see. This inspection interval may be shortened for devices that are moved often or experience unusually heavy use. No other maintenance or service is required.

No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, mechanical, electronic, photocopying, recording, or otherwise, without written permission of Compumedics/Neuroscan.


No license, express or implied, is granted with respect to any of the technology described

in this manual. Compumedics/Neuroscan retains all intellectual property rights associated with the technology described in this manual.

Even though Compumedics/Neuroscan has reviewed this manual, COMPUMEDICS/NEUROSCAN MAKES NO WARRANTY OR REPRESENTATION, EITHER EXPRESS OR IMPLIED, WITH RESPECT TO THIS MANUAL, ITS QUALITY, ACCURACY, MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL COMPUMEDICS/NEUROSCAN BE LIABLE FOR DIRECT, INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES RESULTING FROM ANY DEFECT OR INACCURACY IN THIS MANUAL, even if advised of the possibility of such damage.

## 2 Installation and Details

The MRI Trigger Interface is designed for use with SynAmps<sup>2</sup> systems. It provides an interface between the MR scanner trigger output and the Scan ACQUIRE program. The triggers sent from the scanner will be registered in the ACQUIRE software during acquisition.



### 2.1 Installing the MRI Trigger Interface

The MRI Trigger Interface "splices" in between the Scan computer and the *SynAmps RT* System Unit, using the longer USB cable from the computer to the Interface (USB IN), and the shorter USB cable from the Interface to the System Unit. Alternatively, if you have extra USB connections on the Scan computer, you can connect a USB power cable to the Interface, and use a separate cable to connect the Scan computer to the System Unit.



1. The Interface will accept *optical* or *TTL* trigger pulses, depending on the output generated by your scanner. For example, Siemens has optical trigger outputs, and GE and Varian have TTL pulses. (On GE scanners, use the RF UNBLNK output). *Optical* pulses should be connected to the **MRIIN** input. *TTL* pulses should be connected to the **TDIIN** (Triggering Device Input) input.
2. The MRI Trigger Interface is powered by a USB cable from the *Scan computer*. Connect the USB cable from a USB port on the computer to the **USB IN** input on the Interface.
3. Connect the short USB cable from the **USB OUT** port to the USB port on the back of the *SynAmps RT* System Unit.
4. You can connect the output trigger cable to either the **TriggerOUT** or **Stretched TriggerOUT** connector. **TriggerOUT** passes the trigger directly through, without any changes to the signal (used primarily for debugging). Generally, you will want to extend, or stretch, the duration of the pulse, and it is recommended that you use the **Stretched TriggerOUT** output. (Pulses are stretched to approximately 65 to 70 $\mu$ sec). Connect the other end of the trigger cable to the BNC connection on the *STIM-to-SCAN cable*.
5. TRUE / INVERT switch. In the **TRUE** state, the resting TTL is 5V, and the trigger occurs when the voltage drops to 0V. In the **INVERT** state, the resting voltage is 0V, and trigger occurs at 5V. Set the switch according to your scanner's logic.

## 2.2 In Operation

When connected correctly, you will see the **Pulse** light blink when a trigger is passed through the interface. In ACQUIRE, in the CNT file, you will see a red (response) Type Code of 5 for each pulse from the scanner.

## 2.3 Troubleshooting

If you do not see any trigger events:

1. Verify the type of pulse being sent from the scanner (optical or TTL), and verify it is connected into the correct input on the Interface.
2. Verify all cable connections.
3. If you are not using the Stretched Trigger Output, try using it.
4. Try switching between TRUE and INVERT.

If you are still having problems, contact techsup@neuroscan.com.

If you notice **any** Mechanical Damage, contact techsup@neuroscan.com for return instructions.

DO NOT ATTEMPT TO REPAIR. WARRANTY WILL BE VOID IF UNAUTHORIZED REPAIR WORK HAS BEEN CARRIED OUT.